

Robot Manual (RM_001)

Title: LabAssist Robot Basic Operation

Version: 1.0

Effective Date: 2025-01-01

Purpose:

To guide users in operating and maintaining the LabAssist robot for general laboratory assistance tasks.

Overview:

The LabAssist robot assists with carrying small loads, responding to voice commands, and providing basic information retrieval.

Operation Instructions:

1. Power On: Press the main power button located at the rear until the indicator LED turns green.
2. Connection: The robot connects automatically to the lab Wi-Fi (SSID: LabNet). If connection fails, restart.
3. Commands: Speak clearly using phrases such as "LabAssist, fetch the calibration log" or "LabAssist, summarize SOP_002."
4. Charging: Dock the robot to the charging station when battery level is below 20%.
5. Troubleshooting: Restart if voice commands are not recognized. Contact IT if persistent issues occur.

Notes:

The LabAssist robot is designed to complement human workflows, not replace them.